

URO*logic*

Multi-functional Urological Table





MOVEMENT CAPABILITIES

All the movements of the system are motorized. Tilting movement of the table provides comfort to the user with its -16° (optional, upto -20°) flexibility. Trendelenburg movement allows a range of $+90^\circ$, which is a must for some urological applications.

Placing the patient on the table is even made easier thanks to the wide range of up-down movement. The surface of the table can also move in the transverse and longitudinal direction to avoid moving the patient unnecessarily. Crane movement of the x-ray tube contributes to applying necessary X-ray procedures including radiography without moving the patient.






PCK UROlogic is a fine combination of high technology, versatility and economy. Its ergonomic design allows the user to work without being distracted by any moving parts of the table. Upgradeability provides the user with a variety of choices including the option of adding new features depending on his needs. **PCK UROlogic** determines new standards for urological applications.



BASIC FEATURES

- ▶ High-power X-ray generator
- ▶ High resolution 9" Image Intensifier (12" II, optional)
- ▶ -16° (optional, upto -20°), +90° tilting movement of the table
- ▶ Accurate diagnosis
- ▶ All table accessories needed for all urological applications
- ▶ High quality fluoroscopy and radiography
- ▶ Zoom function
- ▶ Mobile control panel for all X-ray functions
- ▶ High resolution monitor on swiveling arm
- ▶ Emergency stop button
- ▶ Automatic brightness and contrast control
- ▶ Image rotation
- ▶ Foot switch for table movements and fluoroscopy
- ▶ Flexible patient positioning
- ▶ High quality materials
- ▶ Radiation reduction
- ▶ Optimum image quality
- ▶ Ergonomic design for the comfort of user and the patient
- ▶ X-ray system coverage of the whole urinary tract
- ▶ Easy to operate
- ▶ Surface designed for easy clean-up and hygiene
- ▶ The degree of tilting and trendelenburg movements and the height of up-down movements of the table are displayed digitally
- ▶ X-ray warning lamp



ACCESSORIES

- ▶ Accessory stand
- ▶ Standard leg holder
- ▶ Infusion arm board
- ▶ Shoulder supports
- ▶ Elbow supports
- ▶ Grab rails
- ▶ Compression belt
- ▶ Reusable drainage system
- ▶ Disposable bag drainage system
- ▶ Serum hanger
- ▶ Patient covering cloth holder
- ▶ Micturition seat
- ▶ Endoscopy shelf
- ▶ Paper towel hanger
- ▶ Patient support cushion



OPTIONAL FEATURES

- ▶ DSA (Digital Subtraction Angiography)
- ▶ Bucky device
- ▶ Double monitor
- ▶ 12" Image Intensifier
- ▶ Magneto-optical disc archiving system
- ▶ Acquiring, storing and processing of images
- ▶ Sequence recording
- ▶ Archiving facilities
- ▶ Data transfer facilities
- ▶ Patient data base
- ▶ Computer controlled fluoroscopy techniques with real-time digital image processing
- ▶ Multiple patient positioning memories that can be recalled
- ▶ Uro videoscapy

Technical specifications are subject to change without notice.



PCK Electronic Industry and Trade Co., Ltd.

Turan Güneş Blv., Konrad Adenauer Cad. 59/1 Sancak, Çankaya 06550 Ankara TURKIYE

Phone: +90 (312) 491-6010 Fax: +90 (312) 491-6011

PCK Factory Division:

Ankara Sanayi Odası 1. Organize Sanayi Bölgesi Orhan Işık Cad. No: 4 Sincan ANKARA

Phone: +90 (312) 267-2046 Fax: +90 (312) 267-0609

www.pckmed.com

intlsales@pckmed.com

national@pckmed.com